

A New Generation of Suture for Deep Tissue Approximation

The DFX Tensile Tape System



Introducing a new ultra-high molecular weight polyethylene suture for deep tissue approximation.

Imagine a 2mm footprint for greater tissue coverage and hold. Now imagine greater tensile strength, pliability and precise knot placement. Introducing the new DFX Tensile Tape Suture System that leaps beyond conventional deep tissue closure.

The DFX Tensile Tape Suture combines ultra-high molecular weight polyethylene thread with a proprietary flat-braid process. Together, they deliver a low profile, lightweight, non-absorbable platform with high strength and knot pull.

The pliability of the Tape offers excellent passing and coverage when reconstruction of soft-tissue symmetry, deep fascial closure following joint reconstruction and trauma are required.

Discover a new generation of deep tissue approximation. Experience the DFX Tensile Tape System.

DUNAMIS MEDICAL™


877.454.2186
dunamismedical.com




Whether it's an Achilles, Patellar tendon or Quad rupture, the DFX Tensile Tape delivers versatility and confidence.

- The Tape Sutures are offered in 18" and 30" lengths in 2.0mm and 1.5mm widths with two needle options.
- 2.0mm Tensile Tape has a C-7 Half Circle Taper Point Needle that passes beautifully through any deep fascia while approximating muscle, tendons or ligaments.
- 1.5mm Tensile Tape offers a smaller HC-5 Half Circle Needle where precise closure is needed in tight spaces.
- The proprietary flat-braid design offers pliability, increased tensile strength and consistent pressure across the repair site.

Knot Pull Load to Failure (LBF)



Straight Pull Force to Failure (N)



1.5 mm DFX Tensile Tape Suture

Product No.

Box of 12

DFX-0025-S	DFX Tensile Tape 1.5 mm (UHMWPE) 18" White	1/2 Circle Taper Point Needle - 25.9mm HC-5
DFX-0024-S	DFX Tensile Tape 1.5 mm (UHMWPE) 30" White	1/2 Circle Taper Point Needle - 25.9mm HC-5
DFX-0027-S	DFX Tensile Tape 1.5 mm (UHMWPE) 18" Blue	1/2 Circle Taper Point Needle - 25.9mm HC-5
DFX-0026-S	DFX Tensile Tape 1.5 mm (UHMWPE) 30" Blue	1/2 Circle Taper Point Needle - 25.9mm HC-5

2.0 mm DFX Tensile Tape Suture

Product No.

Box of 12

DFX-0021-S	DFX Tensile Tape 2.0 mm (UHMWPE) 18" White	HALF Circle Taper Point Needle - 38.9mm C-7
DFX-0020-S	DFX Tensile Tape 2.0 mm (UHMWPE) 30" White	HALF Circle Taper Point Needle - 38.9mm C-7
DFX-0023-S	DFX Tensile Tape 2.0 mm (UHMWPE) 18" Blue	HALF Circle Taper Point Needle - 38.9mm C-7
DFX-0022-S	DFX Tensile Tape 2.0 mm (UHMWPE) 30" Blue	HALF Circle Taper Point Needle - 38.9mm C-7

Discover a new generation of deep tissue approximation.

877.454.2186
dunamismedical.com

DUNAMIS MEDICAL™