


ORTHOPEDIC

# Discover a better repair.

Rock-solid bite with multiple suture-tape combinations.

## The PUNCHTac™ vPEEK Anchor System



**The PUNCHTac™ 5.5mm** Fully Threaded Anchor introduces an advanced thread pitch design for the shoulder that's a game changer for medial row repair.

It's about maximizing bone-to-anchor surface contact. And it's about distributing stress or torsional load evenly. All of this results in a Pull-Out Force of 554 N<sup>1</sup>.

All anchors are manufactured with Vestakeep® PEEK from Evonik. Vestakeep® PEEK is a biomaterial which exhibits superior modulus and high ductility.



**DUNAMIS** MEDICAL™

Discover a better repair.  
Experience the PUNCHTac™ vPEEK Anchor System.

877.454.2186  
dunamismedical.com




What about footprint coverage or a broader tendon-to-bone contact area?  
The PUNCHTac™ 5.5mm anchors are preloaded with the following tape combinations:

- Double 2.0mm DFX Tensile Tape Sutures, gives you four times standard footprint coverage
- Double #2 Force DFX Suture
- Double #2 Force DFX Sutures and a 1.5mm Tensile Tape

For repairs that present poor bone quality Dunamis Medical offers a triple loaded 6.5mm Fully Threaded Anchor. It's preloaded with two #2 Force DFX Sutures and one 1.5mm DFX Tensile Tape.

### Average Load-to-Failure<sup>1</sup> - Soft Tissue Anchors



### PUNCHTac™ vPEEK Anchor System

Product No.

Box of 5

|          |  |   |
|----------|--|---|
| DFX-0028 | 3.0mm PunchTac™ Anchor vPEEK           | Preloaded with One #2 DFX Force Suture (UHMWPE)                                   |
| DFX-0029 | 3.5mm Salvage PunchTac™ Anchor vPEEK   | Preloaded with Two #2 DFX Force Suture (UHMWPE)                                   |
| DFX-0030 | 4.5mm Knotless PunchTac™ Anchor vPEEK  | Preloaded with One #2 DFX Force Suture (UHMWPE)                                   |
| DFX-0031 | 5.5mm PunchTac™ Threaded Anchor vPEEK  | Preloaded with Two #2 DFX Force Suture (UHMWPE)                                   |
| DFX-0035 | 5.5mm PunchTac™ Threaded Anchor vPEEK  | Preloaded with Two 2mm DFX Tensile Tape (UHMWPE)                                  |
| DFX-0036 | 5.5mm PunchTac™ Threaded Anchor vPEEK  | Preloaded with Two #2 DFX Force Suture<br>and One 1.5mm DFX Tensile Tape (UHMWPE) |
| DFX-0037 | 6.5 mm PunchTac™ Threaded Anchor vPEEK | Preloaded with Two #2 DFX Force Suture<br>and One 1.5mm DFX Tensile Tape (UHMWPE) |

Discover a better repair.  
The PUNCHTac™ vPEEK Anchor System

877.454.2186  
dunamismedical.com

**DUNAMIS** MEDICAL™

<sup>1</sup> Test Data on file at Dunamis Medical  
Vestakeep PEEK is a registered trademark of Evonik®.  
Healicoil PK, TwinFix Ultra HA and TwinFix Ultra PK are trademarks of Smith & Nephew®  
Corkscrew FT II is a trademark of Arthrex®.  
PunchTac is a trademark of Dunamis Medical, LLC  
Listed products are covered by one or more U.S. patents and trademarks.  
Lit#: 03072017P