


ORTHOPEDIC

# Make the switch to better design.

ACL and PCL Reconstruction

## TIGHTFix™ vPEEK Interference Screw System



The TIGHTFix™ vPEEK Interference Screws are a new generation of bone-patella tendon-bone and soft-tissue graft fixation. With an advanced barrel and thread design we've improved bone-to-graft surface contact.

The Insertion Driver tip runs the length of the screw providing axial support. The screws are made with Vestakeep® PEEK, a high strength, bio-inert polymer that exhibits superior modulus and high ductility. The combination of the structural support designed into the driver/screw interface and Vestakeep® that resists fracturing results in greater torsional strength and a stiffer repair.

Also, the TIGHTFix™ ACL-PCL Instrument Set includes a separate tray of taps offering the surgeon intraoperative flexibility.

**DUNAMIS** MEDICAL™


Make the switch to better design.  
Experience the TIGHTFix™ vPEEK Interference Screw System.

877.454.2186  
[dunamismedical.com](http://dunamismedical.com)



TIGHTFix™ vPEEK Interference Screws offer the surgeon an advanced barrel and thread design that improves bone purchase, and balances stress distribution evenly for greater torsional strength.

Average Push-Out Force<sup>1</sup> - Interference Screws



## TIGHTFix™ vPEEK Interference Screw System

Product No.

Each

DFX-0001-T	TIGHTFix™ vPEEK Interference Screw, 8mm x 25mm
DFX-0002-T	TIGHTFix™ vPEEK Interference Screw, 9mm x 25mm
DFX-0003-T	TIGHTFix™ vPEEK Interference Screw, 10mm x 25mm
DFX-0004-T	TIGHTFix™ vPEEK Interference Screw, 11mm x 25mm

Make the switch to better design.  
Experience the TIGHTFix™ System.

877.454.2186  
dunamismedical.com

**DUNAMIS** MEDICAL™

<sup>1</sup> Test Data on file at Dunamis Medical  
Vestakeep® PEEK is a registered trademark of Evonik  
BioScrew is a registered trademark of Linvatec®  
BioRCI and Biosure are trademarks of Smith & Nephew®  
Bio-Interference is a trademark of Arthrex®  
TightFix is a trademark of Dunamis Medical, LLC  
Listed products are covered by one or more U.S. patents and trademarks.  
Lit#: 11266102T