

4.75mm

Double-loaded  
Threaded Anchor



Anchor with Two Fiber-Connect Tapes

Tool Used for  
Implantation


**DFX0043A**

Universal Punch  
DM-4818

4.5/4.75mm Punch Tap



**ANCHOR IMPLANTED  
BY THE SCREW-IN METHOD**



**SUGGESTED USE**

Shoulder-Medial Row, RCR

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)

